

Trinity™ Cup Acetabular System



MODELS

Trinity™ Cup
Trinity™ Polyethylene Liner
Trinity™ ECiMa Liner
Trinity™ Ceramic Liner
Cancellous Bone Screw
Hole Occluder
Trinity™ CoCr Heads
Trinity™ Ceramic Heads

INTENDED PATIENT / USE

The Trinity™ Cup Acetabular System is intended for use in hip replacement for adults (skeletally mature patients). The Trinity™ Cup Acetabular System is intended to provide increased movement and reduce pain by replacing the damaged hip joint.

BENEFITS OF HIP REPLACEMENT

If you have chronic and persistent hip pain that limits regular activity, it is likely you are a good candidate for a hip replacement. As always, it is best to discuss possible treatments with your surgeon.

Pain relief and increased motion:

You should have reduction in pain and better mobility after surgery, within the first few weeks. Many factors, including physical condition, weight, activity level, anatomy and will to comply with your surgeon's instructions before and after surgery will play a key role in your recovery.

PRECAUTIONS

Precautions to take following your hip replacement are dependent on your unique condition and are best advised by your healthcare provider. Seek guidance from a healthcare professional about:

- precautions to take in daily life to guarantee maximum implant survival
- your weight and level of activity can affect the life span of the prosthesis
- you must inform the surgeon of any change in performance (mobility, pain etc.)

POTENTIAL UNDESIRABLE SIDE EFFECTS AND POSSIBLE COMPLICATIONS

Side effects and possible complications related to the placement of the Trinity™ Cup Acetabular System:

- Infection,
- Hip ball and hip socket may separate (hip dislocation),
- Device loosening from the surrounding bone,
- Premature wear or breakage of the implants,
- Bone loss around implant,
- Fracture,
- Hip pain and stiffness,
- Allergic reaction to the implant's materials,
- Change in the length of the treated leg,
- Loss of hip flexibility of the hip joint,
- Audible sounds during motion,
- Nerve damage,
- Damage to blood vessels,
- Embolism,
- Calcification (a process in which calcium builds up causing hardening).

RISK OF DEVICE INTERACTION WITH OTHER EQUIPMENT

The Trinity™ Cup Acetabular System has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. The Trinity™ cup Acetabular System has not been tested for heating, migration, or image artifact in the MR environment.

The safety of the prosthesis in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Typically, airport screening detectors will identify patients with these metal implants. Having a hip implant is common. Let the security officer know that you have a hip replacement and where it is located. Having a metal implant does not prevent you from traveling by air.

EXPECTED DEVICE LIFETIME / INTENDED PERFORMANCE

There are some things that can influence the product lifetime including, but not limited to, surgical indication, surgical technique, patient weight, activity level and comorbidities. You should talk to a registered healthcare professional about your specific situation.

POST OPERATIVE FOLLOW – UP

As determined necessary by a healthcare professional a suitable rehabilitation program should be put in place and regular follow-up checks should be carried out.

Some patients experience swelling of the thigh on the operated side, but this usually disappears quite quickly.

A few patients may experience clicking or other sounds from their new hip, but this rarely causes serious problems and usually disappears after a few months. You must take great care during the first eight to twelve weeks to avoid potentially dislocating your new hip – you must be patient and not try to test your new joint to see how far it will go.

Initially you will tire more easily, not least because there will continue to be traces of anaesthesia in your body for some time. Set aside a rest period each afternoon. You should contact your doctor immediately in the case of any undue pain, severe redness around the operation site or drainage from the wound. Walking without the aid of a stick is often possible from four to six weeks after surgery – although this will be determined by your confidence and progress, and you should follow the advice of your surgeon or physiotherapist.

Your return to driving will be determined by your surgeon which may be as much as six to twelve weeks. Your return to work will also be determined by your surgeon.

Any symptoms that extend beyond the above description should be discussed with a registered healthcare professional. For example: hip pain and stiffness, loss of flexibility of the hip joint, audible sounds during motion, hip dislocation, swelling, fever.

MATERIALS

The Trinity™ Cup, Cancellous Bone Screw and Hole Occluder are manufactured from titanium alloy. The Trinity™ Cup has a two layer coating of calcium phosphate (CaP) on top of titanium.

The Trinity polyethylene liners are manufactured from Ultra High Molecular Weight Polyethylene (UHMWPE) plastic. The Trinity ECiMa liners are manufactured from Ultra High Molecular Weight Polyethylene (UHMWPE) with vitamin E additive.

The Trinity Ceramic Liner and Ceramic Head are manufactured from BIOLOX® delta ceramic.

The Trinity CoCr Head is manufactured from cobalt-chromium alloy.

MANUFACTURING RESIDUALS THAT COULD POSE A RISK TO THE PATIENT

The Trinity™ Cup Acetabular System is cleaned using validated cleaning processes taken from industry standards to ensure that any residuals from the manufacturing process are removed. Testing has demonstrated that there are no residuals that pose a particular risk to the patient.

REPORTING ANY SERIOUS INCIDENT

In the event of a serious incident related to the device contact:

- Your Healthcare Professional
- The manufacturer (Corin Ltd)

Your local Corin representative may also be contacted.

Additionally, your local health authority should be contacted. Refer to the list provided on the back pages of this information leaflet for contact details.

Corin Limited (Manufacturer)

The Corinium Centre
Cirencester,
Gloucestershire
GL7 1YJ,
United Kingdom

Telephone: +44 (0) 1285 659866

Fax: +44 (0) 1285 658960

email: enquiries@coringroup.com

Website: <https://vigilance.coringroup.com/vigilance>

www.coringroup.com/patients/



LOCAL HEALTH AUTHORITY LIST

Australia

The Therapeutic Goods Administration (TGA)

Website: <http://www.tga.gov.au>

België / Belgique / Belgien / Belgium

Federal Agency for Medicines and Health Products (FAMHP)

Website: <http://www.fagg-afmps.be/>

България / Bulgaria

Bulgarian Drug Agency

Website: <https://www.bda.bg/bg/>

Ceska Republika / Czech Republic

State Institute for Drug Control, Medical Devices Branch

Website: <http://www.niszp.cz/>

Danmark / Denmark

Danish Medicines Agency

Website: <https://laegemiddelstyrelsen.dk/da/udstyr/>

Deutschland / Germany

Federal Institute for Drugs and Medical Devices

Website: https://www.bfarm.de/DE/Home/_node.html

Eesti / Estonia

Health Board, Medical Devices Department

Website: <https://www.terviseamet.ee/et>

Eire / Ireland

Health Products Regulatory Authority (HPRA)

Website: <https://www.hpra.ie/>

Ellada / Greece

National Organization for Medicines

Website: <https://www.eof.gr/web/guest/home>

España / Spain

Agencia Española de Medicamentos y Productos Sanitarios

Website: <https://www.aemps.gob.es/>

France

Agence Nationale de Sécurité du Médicament et des Produits de Santé (ANSM)

Website: <https://ansm.sante.fr/>

Hrvatska / Croatia

Agency for Medicinal Products and Medical Devices

Website: <https://www.halmed.hr/?ln=en>

Iceland / Island

Icelandic Medicines Agency

Website: <https://www.lyfjastofnun.is/>

Italia / Italy

Ministry of Health Directorate General of Medical Devices and Pharmaceutical Services

Website: https://www.salute.gov.it/portale/p5_11.jsp

Japan

Pharmaceuticals and Medical Devices Agency (PMDA)

Website: <https://www.pmda.go.jp/english/>

Kypros / Kibris / Cyprus

Cyprus Medical Devices Competent Authority

Website:

<https://www.moh.gov.cy/moh/mphs/mpms.nsf/All/A82FE3D75F4BF2CAC225850A0036075A?OpenDocument>

Latvija / Latvia

Medical Device Evaluation Department State Agency of Medicines

Website: <https://www.zva.gov.lv/en>

Liechtenstein

Office of Public Health

Website: <https://www.llv.li/inhalt/1908/amtsstellen/amt-fur-gesundheit>

Lietuva / Lithuania

The State Health Care Accreditation Agency under the Ministry of Health of the Republic of Lithuania

Website: <https://vaspvt.gov.lt/>

Luxembourg

Ministère de la Santé, Direction de la Santé Division de la pharmacie et des médicaments

Website: <https://sante.public.lu/fr/politique-sante/ministere-sante/index.html>

Malta

Malta Medicines Authority - Medical Devices Unit

Website:

<http://www.medicinesauthority.gov.mt/medicaldevices>

Magyarország / Hungary

National Institute of Pharmacy and Nutrition

Website: https://ogyei.gov.hu/main_page/

Nederland / Netherlands

Ministry of Health, Welfare and Sports CIBG Farmatec-BMC

Website: <https://www.farmatec.nl/>

New Zealand

Medsafe - New Zealand Medicines and Medical Devices Safety Authority

Website - <http://www.medsafe.govt.nz/>

Norway / Norge

Statens legemiddelverk/ Norwegian Medicines Agency

Website: <https://legemiddelverket.no/EnglishÖsterreich/>

Austria

Austrian Federal Office for Safety in Health Care (BASG)

Website: <https://www.basg.gv.at/>



Polska / Poland
Office for Registration of Medicinal Products, Medical Devices
and Biocidal Product
Website: <http://www.urpl.gov.pl/pl>

Portugal
Infarmed - National Authority of Medicines and Health
Products
Website:
<https://www.infarmad.pt/web/infarmad/entidades/dispositivos-medicos>

Romania
National Agency for Medicines and Medical Devices of
Romania
Website: <https://www.anm.ro/>

Slovenija / Slovenia
Agency for Medicinal Products and Medical Devices of the
Republic of Slovenia JAZMP
Website: <https://www.jazmp.si/en/>

Slovenska Republika / Slovakia
State Institute for Drug Control, Medical Devices Section
Website: https://www.sukl.sk/hlavna-stranka/english-version?page_id=256

Suomi / Finland
Finnish Medicines Agency Fimea Medical Devices Unit
Website: <https://www.fimea.fi/>

Sverige / Sweden
Swedish Medical Products Agency 'Läkemedelsverket'
Department of Medical Devices
Website: <https://www.lakemedelsverket.se/sv>

Switzerland
Swissmedic, Swiss Agency of Therapeutic Products
Medical Devices Division
Website:
<https://www.swissmedic.ch/swissmedic/de/home.html>

Turkey
Turkish Medicines and Medical Devices Agency (TITCK)
Ministry of Health
Website: <https://www.titck.gov.tr/>

United Kingdom
Medicines & Healthcare products Regulatory Agency (MHRA)
Website:
<https://www.gov.uk/government/organisations/medicines-and-healthcare-products-regulatory-agency>

United States of America
Food and Drugs Administration (FDA)
Website: <https://www.fda.gov/>

